



Quickly monitor and control your fluid flow and pressure to achieve accurate instrument output and maximized system capabilities with QuickStart™ Sensors from IDEX Health & Science

Now you can easily manage flow and pressure across your entire fluidic system, and save crucial problem-solving time by using our dynamic family of in-line sensors. Our experts have applied decades of knowledge, and over five years of extensive life testing and innovation, to deliver the most advanced transducer technology inside compact, plug-and-play sensors. A fusion of modularity, multiplexing, and intelligent sensing make demanding tasks effortless, giving you

exceptional control over every region of your flow path. Each sensor automatically monitors and provides accurate, real-time data with digital output, allowing you to predict failure, mitigate risk of damage, and optimize your system to maintain maximum performance with ease. Super responsive yet small, our modular sensors are simple to integrate into instruments of any size, and can be installed in arrays to deliver essential data in real time.

PRESSURE SENSORS

Immediately detect blockages and mitigate risk — before important samples are compromised with QuickStart™ Pressure Sensors from IDEX Health & Science. Meticulously engineered and broadly tested for precision sensing, our premium QuickStart Pressure Sensors continually monitor system parameters to provide you the information you need to keep your instrument operating reliably. QuickStart Pressure Sensors connect in-line to your system easily, with minimal solution carryover and bubble trapping.

PRESSURE RANGE
.25 – 14 bar absolute

OPERATING TEMPERATURES
+5 °C – +50 °C

S SMALL FOOTPRINT
1.1 x 1.5 in
Smaller than most other models on the market.

APPLICATIONS
IVD, BIO & POC
Chemically compatible with most reagents.

OUTPUT
DIGITAL



SMALL DEVICES, **HUGE RESULTS**

MODULAR FLOW- • THROUGH DESIGN

Each sensing device has been engineered with an adaptive flow-through design to reduce carryover and prevent bubble trapping with a fully swept sensor flow path. Routine maintenance and replacement is fast too, significantly reducing your instrument downtime.

ARRAY • CAPABILITIES

Built-in intelligence allows you to receive fluidic information across your entire instrument when using multiple sensing devices connected simultaneously.

mmi

HIGH **ACCURACY**

Each sensor is fully calibrated, with verified performance, for high-precision plugand-play applications.

REAGENT **COMPATIBILITY**

LOW POWER

CONSUMPTION

Expertly designed with

proprietary processing

devices deliver a high

level of performance

that maintains lower

power consumption.

with an efficiency

temperatures and

reduces overall

architecture, our sensing

Made from chemically inert materials, our sensors allow you to operate with the majority of reagents used in IVD, BIO and POC applications.

STANDALONE OR • MANIFOLD MOUNT

Sensors are available as a standalone inline component, or can be integrated with a manifold using 4 in-lbs of assembly torque.

SMALL FOOTPRINT

Optimized to make the most of a compact design, our sensors work with any instrument, large or small.

FLOW SENSORS

FLOW RATES -1,000 - 1,000 µl/min

> **OPERATING TEMPERATURES** $+5^{\circ}C - +50^{\circ}C$

> > **SMALL FOOTPRINT** 1.6 x 1.0 x 1.2 in

> > > Smaller than most other models on the market.

APPLICATIONS IVD, BIO & POC

Chemically compatible with most reagents.

OUTPUT DIGITAL

> CONTROL FLOW

Easily manage your flow rates and receive accurate, instantaneous data on system performance — saving critical analysis time — with IDEX Health & Science QuickStart Flow Sensors.

> Extensively designed and thoroughly tested for demanding fluid sensing applications, our superior QuickStart Flow Sensors

identify flow rate variances to keep you informed of system

sensitivities and potential problems

that may require immediate attention. QuickStart Flow Sensors are easy-to-


mount with a quick in-line connection, and are optimized for real-time analysis.

> MITIGATE RISK) MEASURE PERFORMANCE

QuickStart™ Pressure Sensors

PART NUMBERS AND ACCESSORIES		
Part #	Description	
I2C PS200F	200 psi Pressure Sensor Standalone Fitting Option	
I2C PS200M	200 psi Pressure Sensor Manifold Option	
I2C PS200F EVAL	200 psi Pressure Sensor Evaluation Kit	
PSCK-I2C	Pressure Sensor I ² C Connection Kit	

SPECIFICATIONS	
Parameter	Specification
Output Signal	Digital
Operating Voltage	5.0 V
Digital Communication Bus	I ² C
Full Scale Pressure Range	0.25 – 14 bar absolute
Accuracy Below Full Scale	< 1% full scale
Repeatability Error from Zero to Full Scale	1% of measured value or 0.05% of full scale (whichever error is larger)
Pressure Detection Response Time	67 ms
Operating Temperature	+5 °C - +50 °C
Ambient Storage Temperature	-30 °C – +100 °C
Proof Pressure	400 psi
Burst Pressure	800 psi




This graph shows the overlaid output from ten pressure sensors monitoring the same fluidic channel. High reproducibility and the capability for simultaneous reading of multiple sensors make the units extremely valuable for instrumentation applications. In this plot, the sensors are responding to a series of pressure increases over the sensor range.

QuickStart[™] Flow Sensors

PART NUMBERS AND ACCESSORIES		
Part Number	Description	
I2C FS1000F	$1,000~\mu L$ Flow Sensor Standalone Fitting Option	
I2C FS1000M	1,000 µL Flow Sensor Manifold Option	
UART FS1000F	1,000 µL Flow Sensor Standalone Fitting Option	
UART FS1000M	1,000 µL Flow Sensor Manifold Option	
I2C FS1000F EVAL	1,000 µL Flow Sensor Evaluation Kit	
FSCK-I2C	Flow Sensor I ² C Connection Kit	
FSCK-UART	Flow Sensor UART Connection Kit	

SPECIFICATIONS	
Parameter	Specification
Output Signal	Digital
Operating Voltage	3.3 V – 5.0 V
Digital Communication Bus	I ² C / UART
Full Scale Flow Rate	1,000 μL/min
Sensor Output Limit	1,500 µL/min
Accuracy Below Full Scale	5% of measured value or 0.25% of full scale (whichever error is larger)
Repeatability Error from Zero to Full Scale	.5% of measured value or 0.025% of full scale (whichever error is larger)
Flow Detection Response Time	40 ms
Operating Temperature	+5 °C – +50 °C
Ambient Storage Temperature	-40 °C – +80 °C
Proof Pressure	40 psi (3 bar)
Burst Pressure	100 psi (7 bar)



This graph shows the overlaid output from ten flow sensors that are used to monitor the same fluidic channel. In this application, a pump and valve are used to infuse fluid through the fluidic circuit at varying rates. A variety of flow rate changes are observed on the sensors, indicating system response to pump flow and valve switching. Simultaneous monitoring of multiple sensors in such a way can easily show how minute adjustments to system components can affect the measured flow in the fluidic circuit.



For ordering and technical support, please contact:

Europe: Biotech AB Tel +46 (0)300-56 91 80. info@biotech.se www.biotech.se USA: Biotech USA LLC Tel:612-703-5718. sales@biotechusa.us www.biotech.se Japan: BioNik Tel:+81-545-38-9125. info@bionikinc.com www.bionikinc.com