

Analytical and preparative chromatography columns and materials for pharmaceutical applications

○ **Polymeric partition chromatography columns and materials**
MCI GEL™ CHP series

Separation mechanism of CHP series

High performance liquid chromatography relies on one of the following physical phenomena for efficient separation of solutes: partition, adsorption, size exclusion, or ion exchange. Of these, partition chromatography is the most commonly used method, and it separates solutes based on their difference in partitioning between a stationary phase and a mobile phase. This technique has currently become the mainstay in industry for the separation of organic compounds such as pharmaceuticals, agricultural chemicals, and other intermediates. Practically, partition chromatography can be performed in two different modes depending on the relative polarities of the stationary and mobile phases. In the normal phase (NP) mode, the mobile phase is less polar than the stationary phase while the situation is reversed in the reverse phase (RP) mode, where the mobile phase is significantly more polar than the stationary phase.

MCI GEL™ specializes in polymer-based packing materials. The use of polymer-based columns has become more widespread thanks to the many advantages of the polymer matrix like excellent selectivity, the absence of specific adsorption which is found commonly with silica-based packing, operability in a wide pH range and good chemical stability due to the inert nature of polymeric materials. The MCI GEL™ partition chromatography columns are based on a polystyrene and polymethacrylate porous polymer. As RP columns, they are applied to the separation of a wide variety of organic compounds, both in the isocratic and gradient elution mode. The compounds include peptides, insulin, small molecule APIs, nutraceutical compounds, water-soluble vitamins and nucleotides. As NP columns, they are used in the separation of various carotenoids, fat-soluble vitamins, steroids, and food additives. These columns tolerate various organic solvents like hexane, heptane, methylene chloride, and alcohols.

As NP columns, they are used in the separation of various carotenoids, fat-soluble vitamins, steroids, and food additives. Various organic solvents like Hexane Heptane, methylene chloride and alcohols can be used.

The MCI GEL™ packing materials are based on the same chemistries offered in the Diaion™ and Sepabeads™ synthetic adsorbent resins. These polymer chemistries, like Diaion™ HP series and Sepabeads™ SP series, are widely used and documented in the biopharmaceutical industry for fermentation extraction, the food industry and in industrial chromatographic separations. The MCI GEL™ packing materials are available as packed columns for analytical applications, and as bulk packing materials for analytical, preparative and production chromatography applications.

CHP column

MCI GEL™ CHP series are suitable for a wide range of hydrophobicities; porous polystyrene, porous polyacrylate, polymethacrylate. This range of columns is based on the properties of the t

Polystyrene packing: M0

Modified polystyrene pa

Polymethacrylate packing

Modified polymethacryl

The hydrophobicities of the column

MCI GEL™ CHP07/C04, C10 >

Polymer columns for HPLC, with broad pH range, acidic through a

- 1) In reverse phase chromatography, the ionic properties of such compounds would be unsuitable.
- 2) Some extremely hydrophilic compounds like CHP07/C04 or CHP07/C10 could not be separated.
- 3) Polymer columns can be washed with organic solvents.

Polymethacrylate columns, CMC chromatography.


Modified polystyrene packing, interactions occur between the p are difficult to separate using ex phase mode and shows a unique s

All polymeric columns exhibit silanol groups even when end-capped


●CHP column series

Matrix Type	Functionality
Styrene-Divinylbenzene	Non-functional
	Functional

Retentiveness in reverse phase mode



Hydrophobic interaction Interaction mechanism



Durability of polymeric column

The polymeric RP columns are chemically stable. Specifically, the columns have resistance to an alkaline eluent. The following graphs demonstrate stability of the polymeric columns. After feeding a solution of pH 12 into the MCI GEL™ CHP20/C04, there is no change of column performance.

Fig. 5-1 Column durability at pH12 comparison between CHP20/C04 and an ODS column

Conditions
Column : MCI GEL™ CHP20/C04 4.6mm.I.D × 150mm

1,2

1,2

Fig. 5-2 Separation of catecholamines

Conditions
Column : MCI GEL™ CHP20/C04
4.6mm I.D.×150mm
Eluent : 50mM Na-phosphate pH2.5
1.5% Hexanesulfonic acid
CH3CN=80/20
Flow rate : 0.25 ml /min
Column temp. : ambient
Detection : 280nm
Sample : 1. Epinephrine
2. Dopamine
3. 5-Hydroxy tryptophan
4. Serotonin
5. Tryptophan




Fig. 5-4 Purine alkaloids


Conditions
Column : MCI GEL™ CHP20/C04
4.6mm I.D.×150mm
Eluent : H2O/CH3CN=10/90
Flow rate : 0.4 ml /min
Column temp. : 25°C
Detection : 275nm
Sample : 1.Theophylline
2.Theobromine
3.Caffeine

1

2

Application data of CHP series

Fig. 5-6 Uric acid and related compounds



Conditions
 Column : MCI GEL™ CHP20/C04
 4.6mm I.D.×150mm
 Eluent : 10mM TBA/CH₃CN=75/25
 Flow rate : 0.4 mL/min
 Column temp. : 25°C
 Detection : 284nm
 Sample : 1.Hypoxanthine
 2.Xanthine
 3.Theophylline
 4.Uric acid
 5.Orotic acid








Fig. 5-8 Glycyrrhizae ra



Comparison with an ODS column

Fig. 5-7 Bile acids



Conditions
 Column : MCI GEL™ CHP20/C04
 4.6mm I.D.×150mm
 Eluent : 50mM NaH₂PO₄ pH2.0
 /CH₃CN=40/60
 Flow rate : 0.3 mL/min
 Column temp. : 25°C
 Detection : 210nm
 Sample : 1.Cholic acid
 2.Ursodeoxycholic acid
 3.Deoxycholic acid







Fig. 5-9 Adrenal cortex

Application data of CHP series

Fig. 5-10 Ingredients of medicine



Conditions
 Column : MCI GEL™ CMG20/C04
 4.6mm I.D.×150mm
 Eluent : 50mM phosphoric acid(pH2.0)/CH₃OH
 =60/40
 Flow rate : 0.5 mL/min
 Column temp. : 45°C
 Detection : 254nm
 Sample : 1.4-Dimethylaminoantipyrine
 2.Antipyrine
 3.Caffeine
 4.Aspirin
 5.Phenacetin




Fig. 5-12 Peptides

Conditions
 Column : MCI GEL™ CMG20/C04
 4.6mm I.D.×150mm
 Eluent : 0.1%TFA/CH₃CN
 =70/30
 Flow rate : 0.5 mL/min
 Column temp. : 25°C
 Detection : 220nm
 Sample : 1.Gly-Tyr
 2.Met Enkephalin
 3.Leu Enkephalin
 4.Angiotensin II





Fig. 5-14 Procainamide,


Conditions
 Column : MCI GEL™ CMG20/C04
 4.6mm I.D.×150mm
 Eluent : 20mM phosphate pH7
 =65/35
 Flow rate : 0.5 mL/min
 Column temp. : 45°C
 Detection : 254nm
 Sample : 1.Procainamide
 2.Procaine

Comparison with an ODS column

Fig. 5-11 Sulfa drugs



Conditions
 Column : MCI GEL™ CMG20/C04
 4.6mm I.D.×150mm
 Eluent : 20mM phosphate pH6.8/CH₃CN
 =82/18
 Flow rate : 0.5 mL/min
 Column temp. : 45°C
 Detection : 254nm
 Sample : 1.Succinylsulfathiazole
 2.Sulfanilamide
 3.Sulfathiazole
 4.Sulfamerazine



Application data of CHP series

Fig. 5-16 Pravastatin sodium

Conditions
 Column : MCI GEL™ CHP20/C10
 (10μm 250 ×4.6mm I.D.) and
 ODS (10μm 250 ×4.6mm I.D.)
 Eluent : A :0.1% Formic acid;
 B :0.1% Formic acid in AcCN;
 Gradient : 45%B-95%B over 29min.
 Flow rate : 1.00 mL/min
 Column temp.: 25°C
 Detection : UV238nm
 Sample : Pravastatin sodium, Mevastatin and Simvastatin, 1mg/ml each;
 Injection : 5μl





Fig. 5-17 Polyene antibiotics

Conditions
 Column : MCI GEL™ CHP20/C10 (10μm 250 ×4.6mm I.D.) and
 ODS (10μm 250 ×4.6mm I.D.)
 Eluent : A :0.1% Formic acid;
 B :0.1% Formic acid in AcCN; A/B=60/40;
 Flow rate : 1.00 mL/min
 Column temp.: 25°C
 Detection : UV305nm for Nystatin, VIS405nm for
 Amphotericin B and UV340nm for Filipin;
 Sample : Pravastatin sodium, Mevastatin and
 Simvastatin, 1mg/ml each;
 Injection : 10μl





Fig. 5-18 Proteins

Conditions
 Column : MCI GEL™ CMG20/
 4.6mm I.D.×250mm
 Eluent : A 0.05% TFA/CH₃Cl
 B 0.05% TFA/CH₃Cl
 A → B 45min linear g
 Flow rate : 0.5 mL/min
 Column temp.: 25°C
 Detection : 280nm
 Sample : 1. Ribonuclease A
 2. Cytochrome C
 3. Transferrin
 4. α-Chymotrypsin
 5. β-Lactoglobulin

Fig. 5-19 Proteins

CHP20/C10 10μm

Ribonuclease A Co



Application data of CHP series

Fig. 5-20 Insulin

Conditions
 Column : MCI GEL™ CHP20/C10
 MCI GEL™ CMG20/C10
 ODS 10μm
 4.6mm I.D.×150mm
 Eluent : A) 0.1%TFA, H₂O
 B) 0.1%TFA, CH₃OH
 Gradient : 20%B→60%B over 20min.
 Flow rate : 1.0 ml /min
 Column temp.: 40°C
 Detection : 280nm
 Sample : Insulin Glargine and human recombinant , 1mg/ml each
 Injection : 10μl






Fig. 5-21 Ghrelin

Conditions
 Column : MCI GEL™ CMG20/C10
 ODS 10μm
 4.6mm I.D.×150mm
 Eluent : A) 0.1%TFA, H₂O
 B) 0.1%TFA, AcCN
 Gradient : 10%B→60%B over 25min.
 Flow rate : 1.0 ml /min
 Column temp.: 40°C
 Detection : 230nm
 Sample : Ghrelin rat and Ghrelin human ,0.1mmol/l each
 Injection : 10μl

CMG20/C10 10μm

ODS 10μm

Fig. 5-22 Leuprorelin

Conditions
 Column : MCI GEL™ CHP20/C10
 MCI GEL™ CMG20/C10
 ODS 10μm
 4.6mm I.D.×150mm
 Eluent : A) 0.1%TFA, H₂O
 B) 0.1%TFA, AcCN
 Gradient : 20%B→60%B over 20min.
 Flow rate : 1.0 ml /min
 Column temp.: 40°C
 Detection : 280nm
 Sample : Leuprorelin, LHRH human, LHRH salt
 Injection : 10μl




Fig. 5-23 Sifuvirtide

Application data of CHP series

Fig. 5-24 ssRNA Ladder Marker

Conditions
 Column : MCI GEL™ CMG20/C10
 ODS 10μm
 4.6mm I.D.×150mm
 Eluent : A)100mM TEAA, H₂O
 B)100mM TEAA, CH₃CN
 Gradient : CHP10/C10
 ODS 10μm 10%B→40%B over 30min
 Flow rate : 1.0 mL/min 8%B→40%B over 30min
 Column temp.: 40°C
 Detection : 260nm
 Sample : 14-30 ssRNA Ladder Marker [max.0.04mg/ml]
 Injection : 5μl




Fig. 5-25 Nucleotide

Conditions
 Column : MCI GEL™ CHK40/C04
 4.6mm I.D.×150mm
 Eluent : A)19 mM H₃PO₄ / 1 mM NaH₂PO₄ / 5.0% ACN
 B)20 mM Na₂HPO₄ / 100 mM NaClO₄ / 30% ACN
 Gradient : 0-4.0min 0%B 4.0-5.0min 0→30%B 5.0min-6.0min 30%B 6.0min-7.0min 30→50%B
 7.0min-10.0min 50→65%B 10.0min-11.0min 65%B 11.0min- 0%B
 Flow rate : 0.8 mL/min
 Column temp.: 50°C
 Detection : UV260nm
 Sample : 1.Ura, 2.Xan, 3.Thy, 4.Hyp, 5.Gua, 6.Cyt, 7.Ade, 8.Urd, 9.Xao, 10.dT, 11.Ino, 12.Guo, 13.Cyd, 14.Ado
 Injection : 20μl

Urd

Fig. 5-26 Linalool









Fig. 5-27 Coriander

Conditions
 Column : MCI GEL™ CMG20/C10
 4.6mm I.D.×150mm
 Eluent : Hexan/Ethanol=99.1/0.9
 Flow rate : 1.0 mL/min
 Column temp.: 40°C
 Detection : 210nm
 Sample : Coriander
 Injection : 10μl

Application data of CHP series

Fig. 5-28 Application data of CHK40/C04: Separation of Sartans

Conditions
 Column : MCI GEL™ CHK40/C04
 4.6mm I.D.×150mm
 Eluent : A) 10 mM NaH₂PO₄ +0.2 mM Na₂HPO₄ (25%ACN)
 B) 10 mM NaH₂PO₄ +1.0 mM Na₂HPO₄ (40%ACN)
 Gradient : 0.5min 0%B 0.5-2.0min 50%B
 2.0min- 90%B
 Flow rate : 1.0 mL/min
 Column temp.: 50°C
 Detection : UV
 Sample : Candesartan(CAN), Losartan(LOS),
 Valsartan(VAL), Irbesartan(IRB)
 Injection : 20μl



(Data provided by Professor Yokoyama of Yokohama National University)

(Polyphenon 60)

Fig. 5-29 Modified Styrene Divinylbenzene CHP07/C04

Conditions
 Column : MCI GEL™ CHP07/C04
 4.6mm I.D.×150mm
 Eluent : CH₃OH/10mM-Acetic acid=60/40
 Flow rate : 0.46 mL/min
 Column temp.: 60°C
 Detection : 280nm
 Sample : Polyphenon 60(10mg/mL) each 10μL

Fig. 5-30 Styrene Divinylbenzene CHP20/C04

Conditions
 Column : MCI GEL™ CHP20/C04
 4.6mm I.D.×150mm
 Eluent : CH₃OH/10mM-Acetic acid=60/40
 Flow rate : 0.46 mL/min
 Column temp.: 60°C
 Detection : 280nm
 Sample : Polyphenon 60(10mg/mL) each 10μL

(TritonX-100)

Fig. 5-31 C18-alkylated aliphatic

Conditions
 Column : MCI GEL™ CHP07/C04
 4.6mm I.D.×150mm
 Eluent : 50vol%CH₃CN
 Flow rate : 0.50 mL/min
 Column temp.: 40°C
 Detection : 254nm
 Sample : Triton X-100
 (Polyoxyethylene octyl ether)
 1% each 10μL




Fig. 5-33 Application data of

Conditions
 Column : MCI GEL™ CHK45/C05
 4.6mm I.D.×150mm
 Eluent : A) 8 mM H₃PO₄
 B) 10 mM H₃PO₄ /30% ACN
 Gradient : 0-0.7min 0%B 0.7-3.0min 0→40%B 3.0-3.2min 40→80%B 3.2-8.0min 80%B
 Flow rate : 1.3mL/min
 Column temp.: 45°C
 Detection : UV(260nm)
 Injection : 20μl

Analysis of various