

DIAION™ HP20SS

DIAION™ HP20SS is small particle grade based on DIAION™ HP20. A controlled pore size distribution and large surface area offer excellent resolution and the capacity for a wide range of molecules, from small peptides and oligonucleotides up to large proteins. It offers nice balance of pressure flow characteristics and true chromatographic fractionation and has also been successfully applied in simulated moving bed applications for a variety of small bio molecules.

DIAION™ HP20SS is characterized by:

- >> Unique pore size distribution
- >> Excellent batch-to-batch reproducibility
- >> Wide application
- >> High chemical and physical stability
- >> Excellent pressure/flow characteristics

Physical and chemical properties


Grade Name	DIAION™ HP20SS	
Bead form	Spherical, porous	
Matrix	Polystyrene/divinylbenzene	
Chemical Structure	$\begin{array}{c} \text{---CH}_2\text{---CH---CH}_2\text{---CH---} \\ \qquad \qquad \\ \text{C}_6\text{H}_5 \qquad \text{C}_6\text{H}_5 \\ \qquad \qquad \\ \text{---CH---CH}_2\text{---} \end{array}$	
Shipping Density*	g/L	670
Water content	%	55 - 67
Particle Size Distribution on 150 μm	%	15 max.
Particle Size Distribution 63 - 150 μm	%	70 min.
Particle Size Distribution thr. 63 μm	%	20 max.
Particle Density*	g/mL	1.01
Specific Surface Area*	m ² /g	560
Pore Volume*	mL/g	1.2
Pore Radius*	Å	290

Note : properties with a mark "*" are referential data.

Swelling ratio in various solvents

Methanol	1.21
Ethanol	1.21
2-Propanol	1.29
Acetone	1.30
Toluene	1.26
Acetonitrile	1.24
Water	1.00



Pore size distribution**Fig. 1 Pore size distribution of HP20SS****Recommended Operating Conditions**

Maximum Operating Temperature	°C	130
Operating pH Range		0 - 14
Minimum Bed Depth	mm	800
Flow rate	BV/h	Loading 0.5 - 5
	BV/h	Displacement 0.5 - 2
	BV/h	Regeneration 0.5 - 2
	BV/h	Rince 1 - 5
Regenerant		
	Organic solvents for hydrophobic compounds	
	Bases for acidic compounds	
	Acids for basic compounds	
	Buffer solution for pH sensitive compounds	
	Water for an ionic solution	
	Hot steam for volatile compounds	

Hydraulic Characteristics

The approximate pressure drop at various temperatures and flow rates for each meter of bed depth of DIAION™ HP20SS resin in normal down flow operation is shown in the graph below.




Fig. 2 Pressure Drop of HP20SS

FDA status

DIAION™ HP20SS may be used to process food and beverage products and isolate specialized food additives as intended and such used may be said to fully comply with the Federal Food, Drug, and Cosmetic Act.

Applications

- Purification of small peptides, oligonucleotides and proteins
- Adsorption of vitamins, antibiotics, enzymes, steroids and other substance from fermentation solutions
- Decolorization of various sugar solutions
- Adsorption of fatty acids
- Removal of phenol
- Adsorption of various perfume
- Decolorization and purification of various chemicals

Notice

This information are given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. The application, use and processing of our products are beyond our control and therefore your own responsibility.