DIAION™ HP2MGL

DIAION™ HP2MGL is based on crosslinked methacrylate. It does not contain any aromatic compounds. It is considered an intermediate polarity adsorbent resin. It is recommended for desalting and adsorption of organic compounds of relatively high polarity by using the more hydrophilic character of the polymer matrix.

DIAION™ HP2MGL is characterized by:

- >> Unique chemical property and pore size distribution
- >> High chemical and physical stability >> Excellent batch-to-batch reproducibly
- >> Excellent pressure/flow characteristics >> Wide application

Physical and chemical properties

Triyoroan and onemical properties		
Grade Name		DIAION [™] HP2MGL
Bead form		Spherical, porous
Matrix		Crosslinked polymethacrylate
Chemical Structure		$\begin{array}{cccccccccccccccccccccccccccccccccccc$
Whole beads count	-	95 min.
Shipping Density*	g/L	720
Water content	%	55 - 65
Particle Size Distribution thr. 355 μm	%	1 max.
Effective size	mm	0.40 min.
Uniformity Coefficient	-	1.6 max.
Particle Density*	g/mL	1.09
Specific Surface Area*	m²/g	570
Pore Volume*	mL/g	1.3
Pore Radius*	Å	240

Note: properties with a mark "*" are referential data.

Swelling ratio in various solvents

Methanol	1.02
Ethanol	1.05
2-Propanol	1.02
Acetone	1.04
Toluene	1.07
Acetonitrile	1.01
Water	1.00



Pore size distribution




Fig. 1 Pore size distribution of HP2MGL

Recommended Operating Conditions

Maximum Operating Temperature	°C	130		
Operating pH Range		0 - 14		
Minimum Bed Depth	mm	800		
Flow rate	BV/h	Loading 0.5 - 5		
	BV/h	Displacement 0.5 - 2		
	BV/h	Regeneration 0.5 - 2		
	BV/h	Rince 1 - 5		
Regenerant				
Organic solvents for hydrophobic compounds				
Bases for acidic compounds				
Acids for basic compounds				
Buffer solution for pH sensitive compounds				
Water for an ionic solution				
Hot steam for volatile compounds				



Hydraulic Characteristics

The approximate pressure drop at various temperatures and flow rates for each meter of bed depth of DIAIONTM HP2MGL resin in normal down flow operation is shown in the graph below.




Fig. 2 Pressure Drop of HP2MGL

FDA status

DIAIONTM HP2MGL may be used to process food and beverage products and isolate specialized food additives as intended. Such use may be said to fully comply with the Federal Food, Drug, and Cosmetic Act, and applicable food additive regulations, including 21 CFR 177.2470 (Polyester resins, cross-linked).

Applications

- Purification of small peptides, oligonucleotides and proteins
- Adsorption of vitamins, antibiotics, enzymes, steroids and other substance from fermentation solutions
- Decolorization of various sugar solutions
- Adsorption of fatty acids
- Adsorption of various perfume
- Decolorization and purification of various chamicals

Notice

This information are given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. The application, use and processing of our products are beyond our control and therefore your own responsibility.

