

## Product Data Sheet


## SEPABEADS™ SP207

SEPABEADS™ SP207 is modified highly porous styrenic adsorbents. It has higher hydrophobicity and greater selectivity for non-polar molecules, which is derived from chemically bonded bromine to the aromatic rings, than standard aromatic adsorbents. It is suitable for upward flow and batch processes due to its high particle density.

SEPABEADS™ SP207 is characterized by:

- >> Unique chemical structure and higher hydrophobicity
- >> High particle density
- >> Excellent batch-to-batch reproducibly
- >> High chemical and physical stability
- >> Wide application

## Physical and chemical properties

Grade Name	SEPABEADS™ SP207	
Bead form	Spherical, porous	
Matrix	Modified polystyrene/divinylbenzene	
Chemical Structure		
Whole Bead Count	-	95 min.
Shipping Density*	g/L	790
Water content	%	43 - 53
Particle Size Distribution thr. 250 µm	%	10 max.
Effective size	mm	0.25 min.
Uniformity Coefficient	-	1.6 max.
Particle Density*	g/mL	1.18
Specific Surface Area*	m <sup>2</sup> /g	600
Pore Volume*	mL/g	1.0
Pore Radius*	Å	110


Note : properties with a mark "\*" are referential data.

## Swelling ratio in various solvents

Methanol	1.10
Ethanol	1.11
2-Propanol	1.12
Acetone	1.13
Toluene	1.13
Acetonitrile	1.12
Water	1.00



Phone: 212-204-0075  
 Email: [info@pyvot.tech](mailto:info@pyvot.tech)  
 Web: [www.pyvot.tech](http://www.pyvot.tech)

**Pore size distribution****Fig. 1 Pore size distribution of SP207****Recommended Operating Conditions**


Maximum Operating Temperature	°C	130
Operating pH Range		0 - 14
Minimum Bed Depth	mm	800
Flow rate	BV/h	Loading 0.5 - 5
	BV/h	Displacement 0.5 - 2
	BV/h	Regeneration 0.5 - 2
	BV/h	Rince 1 - 5
Regenerant		
	Organic solvents for hydrophobic compounds	
	Bases for acidic compounds	
	Acids for basic compounds	
	Buffer solution for pH sensitive compounds	
	Water for an ionic solution	
	Hot steam for volatile compounds	



# SEPABEADS™ SP207

## Hydraulic Characteristics

The approximate pressure drop at various temperatures and flow rates for each meter of bed depth of SEPABEADS™ SP207 resin in normal down flow operation is shown in the graph below.



**Fig. 2 Pressure Drop of SP207**

## Applications

- Purification of small peptides, oligonucleotides and proteins
- Adsorption of vitamins, antibiotics, enzymes, steroids and other substance from fermentation solutions
- Decolorization of various sugar solutions
- Adsorption of fatty acids
- Removal of phenol
- Adsorption of various perfume
- Decolorization and purification of various chemicals

## Notice

This information are given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. The application, use and processing of our products are beyond our control and therefore your own responsibility.



Phone: 212-204-0075  
Email: [info@pyvot.tech](mailto:info@pyvot.tech)  
Web: [www.pyvot.tech](http://www.pyvot.tech)