

Product Data Sheet


SEPABEADS™ SP207SS

SEPABEADS™ SP207SS is a small size grade based on SEPABEADS™ SP207. It has higher hydrophobicity and greater selectivity for non-polar molecules, which is derived from chemically bonded bromine to the aromatic rings, than standard aromatic adsorbents. It is applied to reversed phase chromatography.

SEPABEADS™ SP207SS is characterized by:

- >> Unique chemical structure and higher hydrophobicity
- >> High performance for chromatography
- >> Excellent batch-to-batch reproducibly
- >> High chemical and physical stability
- >> Wide application

Physical and chemical properties

Grade Name	SEPABEADS™ SP207SS	
Bead Form	Spherical, porous	
Matrix	Modified polystyrene/divinylbenzene	
Chemical Structure		
Shipping Density*	g/L	790
Water Content	%	43 - 53
Particle Size Distribution on 150 µm	%	15 max.
Particle Size Distribution 63 - 150 µm	%	70 min.
Particle Size Distribution thr. 63 µm	%	20 max.
Particle Density*	g/mL	1.18
Specific Surface Area*	m ² /g	590
Pore Volume*	mL/g	1.0
Pore Radius*	Å	110


Note : properties with a mark "*" are referential data.

Swelling ratio in various solvents

Methanol	1.11
Ethanol	1.17
2-Propanol	1.19
Acetone	1.20
Toluene	1.19
Acetonitrile	1.20
Water	1.00



Phone: 212-204-0075
 Email: info@pyvot.tech
 Web: www.pyvot.tech

Pore size distribution**Fig. 1 Pore size distribution of SP207SS****Recommended Operating Conditions**

Maximum Operating Temperature	°C	130
Operating pH Range		0 - 14
Minimum Bed Depth	mm	800
Flow rate	BV/h	Loading 0.5 - 5
	BV/h	Displacement 0.5 - 2
	BV/h	Regeneration 0.5 - 2
	BV/h	Rince 1 - 5
Regenerant		
Organic solvents for hydrophobic compounds		
Bases for acidic compounds		
Acids for basic compounds		
Buffer solution for pH sensitive compounds		
Water for an ionic solution		
Hot steam for volatile compounds		



SEPABEADS™ SP207SS

Hydraulic Characteristics

The approximate pressure drop at various temperatures and flow rates for each meter of bed depth of SEPABEADS™ SP207SS resin in normal down flow operation is shown in the graph below.




Fig. 2 Pressure Drop of SP207SS

Applications

- Purification of small peptides, oligonucleotides and proteins
- Adsorption of vitamins, antibiotics, enzymes, steroids and other substance from fermentation solutions
- Decolorization of various sugar solutions
- Adsorption of fatty acids
- Removal of phenol
- Adsorption of various perfume
- Decolorization and purification of various chemicals

Notice

This information are given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. The application, use and processing of our products are beyond our control and therefore your own responsibility.



Phone: 212-204-0075
 Email: info@pyvot.tech
 Web: www.pyvot.tech