

# Evaluation of Porous Layer Thickness of Core Shell Particle for Separation of Proteins

Norikazu Nagae (1), Tomoyasu Tsukamoto (1), Shun Kojima (1), Scott Silver (2)

1) ChromaNik Technologies Inc., 6-3-1 Namiyoke, Minato-ku, Osaka, Japan


2) Innovations United, 300 East 57th Street Suite 11J, New York, NY 10022, USA

## Abstract




The feature of superficially porous (core shell) particle used as a highly efficient material is existence of a core, a thin porous layer and narrow particle size distribution, which lead to higher efficiency than totally porous particle. Recently a core shell particle with wide pore for biomacromolecular separations has developed by a few manufacturers. It has been said that thin porous layer of core shell particle have an advantage for separation of large molecules such proteins because a diffusion coefficient becomes small to proportional to a molecular weight and a mass transfer speed also decreases. In this paper, thickness of porous layer of core shell particle was evaluated to separate proteins. 2 kinds of thickness of porous layer such as 0.2  $\mu\text{m}$  and 0.5  $\mu\text{m}$  thickness were applied for separation of standard protein samples. On fast separation, 0.2  $\mu\text{m}$  of porous layer showed sharper peaks than 0.5  $\mu\text{m}$  of porous layer. However at 80 degree Celsius and using 60 min gradient time program, 0.5  $\mu\text{m}$  of porous layer showed much sharper peaks than 0.2  $\mu\text{m}$  of porous layer. It was considered that 0.5  $\mu\text{m}$  of porous layer had a wider specific surface area than 0.2  $\mu\text{m}$  of porous layer and this wider specific surface area leaded separation efficiency concerning the partition interaction on the stationary phase to be large.





## Effect of temperature



## Comparison of thickness of porous layer




## Effect of gradient time



Column:  
SunShell C8-30HT, 3.4  $\mu\text{m}$  (30 nm, 0.2  $\mu\text{m}$  layer) 100 or 150 x 2.1 mm  
SunShell C8-30, 2.6  $\mu\text{m}$  (30 nm, 0.5  $\mu\text{m}$  layer) 100 or 150 x 2.1 mm  
Mobile phase: A) 0.1% TFA in water, B) 0.08 % TFA in Acetonitrile  
Gradient program: Time 0 min 5 or 60 min  
%B 20% 65%  
Flow rate: 0.5 mL/min  
Temperature: 80  $^{\circ}\text{C}$   
Detection: UV@215 nm  
Sample: 1 = Cytochrome C, 2 = Lysozyme, 3 = BSA, 4 = Myoglobin, 5 = Ovalbumin

✓ In case of fast separation using 5 minute gradient program,



✓ In case of high resolution separation using 60 minute gradient program at 80  $^{\circ}\text{C}$ ,

reversely is better than