RAPID CRS Evidence Sheet —Real-time PCR—

[Abstract]

Few plate seal for PCR has DNA free and DNase/RNase certificate. Therefore, we developed a certified seal, RAPID CRS. The compatibility of the RAPID CRS with qRT-PCR (SYBR Green Method) was evaluated through comparison of qRT-PCR reactions between two sealing films, RAPID CRS and a seal recommended by a PCR instrument manufacturer. The same result between them showed that the RAPID CRS is useful for qRT-PCR.


[Methods]

Firstly, we generated cDNA from mRNA extracted human normal cells and ran qRT-PCR reactions with each pair of two genes primers. Secondly, we used KAPA SYBR FAST qPCR Kit (Kapa Biosystems, Inc.) and pipetted samples into a 96-well plate. Thirdly, we sealed the right-half side of the plate with RAPID CRS and the left-half side of plate with a seal recommended by a PCR instrument manufacturer. Finally, we ran qRT-PCR by using TaKaRa PCR Thermal Cycler Dice and compared results between them in triplicate.


[Results]

The results showed similar PCR reactions for gene 1 and 2 despite different sealing films (Figure 1, Table 1). Little difference of cycle threshold values in triplicate showed that RAPID CRS is as good as the recommended sealing film.


<Figure 1> Amplification Plots & Dissociation Curve




Manufacturer's Recommended Sealing Film



Comparison of merged data of two results above





Gene1

Samples	Ct (SDM)	Average
RAPID CRS 1	21.29	
RAPID CRS 2	21.35	21.30
RAPID CRS 3	21.25	
Recommended Seal 1	21.32	
Recommended Seal 2	21.32	21.30
Recommended Seal 3	21.26	

Gene2

Samples	Ct (SDM)	Average
RAPID CRS 1	24.81	
RAPID CRS 2	24.75	24.87
RAPID CRS 3	24.78	
Recommended Seal 1	24.82	
Recommended Seal 2	24.86	24.78
Recommended Seal 3	24.92	

SYBR is a registered Trademark of Molecular Probes, Inc.

