

Thank you for purchasing our RAPID SlitSeal.
Please read this manual carefully before use and
save this manual for your future use.

Please note that product specifications and packings may
change without prior notice.

Announcement for your use with ease.

【Purpose of use】

This product can be used for research purposes only.

【Notes for use】

Handle this product carefully when unpacking from the bag.


Do not process this product in autoclave treatment.

【Storage and disposal of product】


Store at room temperature and in a place not exposed to direct
sunlight. To dispose the product after use, comply with local laws
and ordinances, and process properly according to the regulations
of your facility

Instructions


*Please wear gloves before use.



Remove the release liner with a red indication and expose the adhesive
surface.



Remove the protective film before sampling or dispensing. If used without removing the protective film, the needle or the tip
can not puncture the film. If the protective film stays on and is left at 100°C, the edge of the protective film may be curved.



Paste the RAPID SlitSeal on the 96-well plate and make sure
holes of RAPID SlitSeal are correctly positioned on each
well spot.



1040 1st Avenue
Suite # 330
New York, NY 10022
Phone: 212-204-0075
Email: info@pyvot.tech
Web: www.pyvot.tech

| | |
|--------------------------------|-----------------------------|
| Catalog Code | IURSS-S96-801- 22 PV03-1002 |
| Product name | RAPID Slit Seal |
| Bonding method | Pressure sensitive adhesive |
| Materials | PET, Silicone |
| Operating temperature limit | -80°C~100°C * |
| Size | W 80 mm × L 122 mm |
| Quantity | 100 sheets/Box |

* The upper limit of the operating temperature is specified as 100 ° C, and this is clarified by the evaluation
test carried at 100 ° C for 30 minutes.