Evaluation of Durability of Hybrid Silica with Ethylene Chains

N. Nagae^a, T. Tsukamoto^a, R. Koyama^a, E. Shearer^b BioNik Inc. [a] ChromaNik Technologies Inc. Namiyoke, Minato-ku, Osaka Japan (b] BioNik Inc. Obuchi, Fuji, Shizuoka Japan



Silica-based reversed-phase columns have been widely used since the 1970s. Silica has the advantage of being mechanically strong and able to withstand high pressures, but conventional silica C18 has limited use under alkaline conditions and is only durable up to pH 8. In the 2000s, hybrid silica with ethylene chains bonded to the inside or surface of silica was developed, and patent applications were filed for these technologies. Currently, these technologies are freely available, and it is expected that these hybrid silica packings will be available on the market.

In this study, we evaluated the durability of so-called hybrid silica with ethylene chains bonded to the inside and surface of silica, and hybrid C18 bonded with C18 and end-capped. The amount of hybrid silica dissolved was measured when an alkaline mobile phase of pH 11.5 was passed through the column, and the hybrid C18 bonded with C18 and end-capped was also evaluated in the same way. The durability of hybrid silica under alkaline conditions greatly differed depending on the method of introducing and bonding ethylene chains to silica, with a difference of 100 to 1000 times the amount of dissolved silica. The difference in stability of hybrid C18 under alkaline conditions was smaller

than that of hybrid silica alone, suggesting that the effect of end-capping on durability is greater than that of hybrid silica.

2 types of hybrid silica gel **Comparison of durability of hybrid silica gel under alkaline conditions** 2. **Hybrid (2)** 1. Hybrid (1)**Durability test conditions (pH 11.5)** Ethylene chains bonded to the silica gel surface Ethylene cross-linked silica gel pH 11.5, 40 °C, 17 hour Column dimension: 100 x 2.1 mm Mobile phase: 50 mM potassium side, 20 (20 °C) 20 phosphate pH 11.5 Flow rate: 0.1 mL/min (1 mol) (4 mol 16.5 Temperature: 40 °C 15.4 Elution time: 17 h 15 Ð (The carbon loading of the ethylene chain is shown in parentheses.) 10 **Hybrid** (6.4%) The amount of denting increases Full hybrid Company A (7%) exponentian Company B (7.3%) ? Full or partial hybrid of **Hybrid (2**.3%) 1.2 12.5 nm, 125 Å 1.5 Partial hybrid 0.3 (64 h) **Hybrid (**1.2%**)** 6.6 nm, 66 Å Thickness of silica skeleton (41 h) 0.1 Thickness of silica skeleton Silica with only Silica with Hybrid ① No hybrid (340 m²/ end-cappin Hybrid (2) end-capping (4.4%) Fully hybrid silica (4.4%) Partially hybrid silica 3.5 %C: 7% pH 11.5, 17 hour, 40 °C, 60 °C **Durability test conditions (pH 11.5)** %C: 2.5% Column dimension: 100 x 2.1 mm e vithin a few hours, the $(340 \text{ m}^2/\text{g})$ Mobile phase: 50 mM potassium phosphate pH 11.5 olumn Þressure ncreased so much that Flow rate: 0.1 mL/min the measurements were Temperature: 40 °C and 60 °C 2.3 stopped.

(180 m²/g)

Hybrid silica gels used in this study

Name	Specific surface area	Pore diameter	Carbon loading	reagent	
Hybrid ① (6.4%)	180 m²/g	15 nm	6.4%	$\begin{array}{c} OC_{2}H_{5} \\ I \\ H_{5}C_{2}O \\ OC_{2}H_{5} \\ OC_{2}H_{5} \end{array} + \begin{array}{c} H_{5}C_{2}O \\ H_{5}C_{2}O \\ H_{5}C_{2}O \\ H_{5}C_{2}O \end{array} OC_{2}H_{5} \\ OC_{2}H_{5} \\ OC_{2}H_{5} \end{array} + \begin{array}{c} H_{5}C_{2}O \\ H_{5}C_{2}O \\ H_{5}C_{2}O \\ OC_{2}H_{5} \\ OC_{2}H_{5} \end{array} OC_{2}H_{5} \\ OC_{2$	
Hybrid ② (2.3%)	340 m²/g	12 nm	2.3%	Silica gel + H_{3C}	
Hybrid ② (1.2%)	340 m²/g	12 nm	1.2%	Silica gel + $H_{3C} H_{3C} H$	
Silica with end- capping (4.4%)	340 m²/g	12 nm	4.4%	Silica gel + H_3C H_3C H_3C H_3C H_3C H	



- When comparing the amount (speed) of silica dissolution when an alkaline mobile phase was passed through hybrid silica, Hybrid 2 showed a value 100 to 1000 times higher than Hybrid 1.
 The amount of silica C18 dissolved in an alkaline mobile phase increases exponentially over time, and this was also the case for hybrid silica (Hybrid 1).
 When the temperature difference was 20°C, the amount of silica dissolved was about 10 times
- different, and a 10°C temperature increase caused the amount of silica dissolved to increase by about three times.

C18 bonding and end-capping



Durability under basic pH condition Durable test condition Column size: 150 x 4.6 mm Mobile phase: Methanol/50mM potassium phosphate pH 11.5=10/90

Conclusions

 It was confirmed that the durability of hybrid silica under alkaline conditions varies greatly depending on the method of bonding the ethylene chains.

Elution time: 17 h



 It was suggested that the durability of the C18 packing after C18 bonding and end-capping is greatly affected not only by the durability of the hybrid silica, but also by the surface treatment, especially end-capping.

C18 using hybrid silica with ethylene chains bonded within the silica skeleton was found to have sufficient durability under mobile phase conditions from pH 1 to pH 12.

Evaluation of Bidentate End-capping Silvlation Reagents for HPLC

N. Nagae^a, T. Tsukamoto^a, R. Koyama^a, E. Shearer^b [a] ChromaNik Technologies Inc. Namiyoke, Minato-ku, Osaka Japan [b] BioNik Inc. Obuchi, Fuji, Shizuoka Japan



www.chromanik.co.jp



★Six types of C18 packing materials were prepared and compared for hydrogen bonding, hydrophobicity, steric selectivity, peak shape for basic compounds, and durability against acidity and alkalinity.

Tanaka Index, specific surface area and

carbon loading

Opyvol

Phone: 212-204-0075

BioNik Inc.

www.bionikinc.com



Peak shape of basic compounds







- Column: Shown in figure
- Column size: 150 x 4.6 mm Mobile phase: Acetonitrile/20mM phospha
- Mobile phase: Acetonitrile/20mM phosphate buffer pH7.0=60/40 Flow rate: 1.0 mL/min
- Temperature: 40 °C
- Detection: UV@250 nm
- Sample: 1=Uracil, 2=Propranolol, 3= Nortriptyline, 4=Amitriptyline



N: Theoretical plate, TF: Tailing factor

★Under neutral conditions using acetonitrile as the mobile phase, basic compounds are prone to tailing. With bidentate end-capping, both types A and B showed less tailing. The results were more symmetrical than Company W Hybrid C18, which was used for comparison.

Durability under acidic and basic pH conditions



Durable test condition Column size: 50 x 2.1 mm Mobile phase: Acetonitrile/1.0% TFA, (pH 1)=10/90 Flow rate: 0.4 mL/min Temperature: 80 °C

Measurement condition Column size: 50 x 2.1 mm

Conclusions

- Six types of packing materials were synthesized using three types of silica base materials and two types of bidentate end-capping reagents after C18 bonding.
- Comparison of hydrogen bonding, hydrophobicity, and steric selectivity in columns packed with six types of packing materials showed that the endcapping reagent had little effect, while the silica base material had a greater effect.

Mobile phase: Acetonitrile/water =60/40 Flow rate: 0.4 mL/min, Temperature: 40 °C Sample: 1 = Uracil (t₀), 2 = Butylbenzene

Durable test condition Column size: 150 x 4.6 mm Mobile phase: Methanol/50mM potassium phosphate pH 11.5=10/90 Flow rate: 1.0 mL/min Temperature: 40 °C

Measurement condition Mobile phase: Acetonitrile/water=70/30 Flow rate: 1.0 mL/min Temperature: 40 °C Sample: 1 = Butylbenzene

★ For comparison, the test results of hybrid C18 from companies W, Y, and P are also shown.

The degree of tailing of the peaks of basic compounds was compared. As a result, no significant difference was observed among the six types of packing materials, and peaks with less tailing were obtained. Company W hybrid C18, used as a reference, showed tailing of basic compounds (especially nortriptyline).

✓ In both durability tests under acidic and basic pH conditions, packing materials using B type 1,2-bis(chlorodimethylsilyl)ethane for end-capping were more durable than A type 1,5-dichlorohexamethyltrisiloxane. It was presumed that this is because the B type end-capping reagent is more hydrophobic, making it difficult for water molecules, which cause hydrolysis and lead to deterioration, to approach the packing silica surface.

The durability under basic pH condition was greatly affected by the silica base material, and the silica with ethylene chains crosslinked into the silica skeleton showed the highest durability, which was equivalent to that of Hybrid C18 from company W.

Evaluation of Hybrid Silica C18 End-capped with Bidentate Silylating Reagent for HPLC Opyvot Phone: 212-204-0075 Email: info@pyvot.tech

N. Nagae^a, T. Tsukamoto^a, R. Koyama^a, E. Shearer^b Bio Nik Inc. [a] ChromaNik Technologies Inc. Namiyoke, Minato-ku, Osaka Japan [b] BioNik Inc. Obuchi, Fuji, Shizuoka Japan www.bionikinc.com

ChromaNik Technologies Inc.

www.chromanik.co.jp

Comparison between Company W hybrid C18 and SunBridge C18 with Ethylene cross-linked silica gel And 1,2-Bis(chlorodimethylsilyl)ethane end-capping









Conclusions



- ✓ SunBrige C18 was developed using Ethylene cross-linked silica gel and 1,2-Bis(chlorodimethylsilyl)ethane end-capping reagent.
- ✓ As a comparison, Hybrid C18 manufactured by W Company, which is ethylene cross-linked silica gel C18 in which an ethylene chain is incorporated into the silica skeleton, was used.
- \checkmark Company W Hybrid C18 showed severe tailing of metal-chelating, oxine. Compared to Company W Hybrid C18, SunBridge C18 showed a higher number of plates and less tailing of peaks for metal-chelating and basic compounds.
- ✓ SunBridge C18 showed almost the same or higher stability under pH 11.5 and pH 1.0 conditions as Company W Hybrid C18, which shows an order of magnitude higher stability than other companies' hybrid C18.
- ✓ SunBridge C18 showed a excellent reproducibility in retention factor under 100% aqueous conditions while Company W Hybrid C18 showed a poor reproducibility.

Analysis of PFAS in Tap Water Using a Pentafluorophenyl Column Oyvo



Tadashi Kitta¹, Hirotake Takahashi¹, Norikazu Nagae², Ryuji Koyama² and Tomoyasu Tsukamoto² ¹Japan Food Inspection Corporation, Heiwajima, Ohta-ku, Tokyo, Japan ²ChromaNik Technologies Inc., Namiyoke, Minato-ku, Osaka, Japan

ChromaNik Technologies Inc.

Perfluoroalkyl and polyfluoroalkyl substances (PFAS) are known to pose health risks due to their bioaccumulative nature and environmental persistence. In April 2026, legal restrictions setting a combined limit of 50 ng/L for PFOS and PFOA in tap and mineral water will take effect in Japan. Furthermore, the European Union has already established maximum levels for four PFAS congeners in seafood, and similar regulations are anticipated in the United States.

In this study, we developed a pentafluorophenyl (PFP) column consisting of an ethylene cross-linked hybrid silica gel and a PFP stationary phase. This column is capable of simultaneously analyzing 19 PFAS compounds, ranging from short-chain to long-chain substances. Its applicability to drinking water analysis was evaluated and compared with that of a conventional C18 column. The PFP column exhibited not only greater sensitivity for PFBA but also superior separation of PFHxS isomers compared to the C18 column.

Hybrid silica PFP and C18

Per- and polyfluoroalkyl substances (PFAS) are defined as "fluorinated substances that contain at least one fully fluorinated methyl or methylene carbon (with no H/Cl/Br/l atoms attached). As of January 2025, there are 14,735 PFAS

with clear chemical structures and 1,915 PFAS with unclear structural information such as polymers and complex mixtures.



PFAS are heat-resistant, water-repellent, and oil-repellent, and are used in many industrial products and daily necessities. On the other hand, they are called "Forever Chemicals" because they are difficult to decompose in the environment and higher bioaccumulative. The main routes of human exposure are oral ingestion of drinking water and food, inhalation of product dust, and skin absorption through contact with products (such as cosmetics).



Analysis issues When a conventional C18 column is used, 1) lack sensitivity for short-chain PFAS. 2) insufficient separation of long-chain PFAS isomers. \rightarrow A method that can comprehensively and sensitively analyze a variety of PFASs is needed.

Objectives Aim to establish a comprehensive and highly sensitive method for analyzing trace PFAS in tap water using PFP columns.

Flow rate

(mL/min)

0.30

0.30

0.30

0.40

0.40

0.40

0.40

0.30

%B

30

30

90

97

97

30

30

30

Measurement conditions	

HPLC	Ultimate 3000(Thermo Fisher Scientific)

Sample pretreatment

Comparison of isomer separation

NL: 1.20E6





Comparison of PFBA peak intensity between PFP column and C18 column

Column	Peak intensity				
	100% MeOH	30% MeOH +70% ammonium formate	100% MeCN		
SunBridge PFP-R	442,000	173,000	336,000		
SunBridge C18	277,000	110,000	218,000		

Isocratic separation using 100% MeOH, 100% MeCN and 30% MeOH + 70% ammonium formate. Sample: PFBA (5 ng/mL), injection volume: 4 μ L



Comparison of separation of 10 kinds of PFASs







Conclusions

- ✓ SunBridge PFP-R exhibited not only greater sensitivity for PFBA but also superior separation of PFHxS isomers compared to SunBridge C18.
- ✓ SunBridge PFP-R was able to completely separate 19 kinds of PFASs although the C18 was unable to adequately separate PFNA and PFOS.
- ✓ SunBridge PFP-R is made of hybrid silica gel, so it is expected to have a longer column life than conventional PFP columns.