

DGC-020-153881

Product Data Sheet

Product number: 115388100001

FULLY SWEEP
INLINE
PATHWAY

TROUBLE
FREE FLUID
CONNECTION

ACTIVE
DISSOLVED GAS
REMOVAL

High-Efficiency Degassing for Your Fluidic System

Improve data quality and reliability in your system with degassing chambers from IDEX Health & Science.

IDEX Health & Science inline degassing chambers provide superior fluid conditioning for life science instrumentation. Our degassing chambers improve instrument precision and assay quality by removing dissolved gases from aqueous system fluids before they outgas and form bubbles that impact precision dispensing, prevent mixing, disrupt separation air gaps, or interfere with detection.

Features:

- › Active dissolved gas removal
- › Fully swept inline fluid path
- › Trouble free fluid connection with 1/4-28 UNF fittings or manifold mountable
- › Vacuum connection for elastomeric tubing
- › Sturdy housing for membrane protection

Applications:

- › Benchtop clinical chemistry
- › Benchtop immunoassays
- › Flow cytometry
- › NGS applications
- › Spatial biology
- › Proteomics

Description

The core functional element of this inline degasser is a fluid path lined by a highly permeable silicone-like membrane. Fluid flows through the membrane pathway, while the space surrounding the membrane is evacuated. A PID-controlled IDEX Health & Science vacuum pump holds a precise vacuum level driving the active removal of dissolved gases across the membrane.



Phone: 212-204-0075
Email: info@pyvot.tech
Web: www.pyvot.tech



Intelligent Solutions for Life®
Fluidics | Optics | Consumables | Assemblies

General

Parameter	Value
Class	Degassing chamber
Series	DGC
Product alias	DGC-020-153881
Product number	115388100001

Absolute Maximum Ratings (per channel)

Parameter	Value	Unit
Maximum operating pressure difference between fluid and vacuum	100	kPa
Maximum operating flow rate	20 333	mL/min µL/sec
Maximum operating temperature	40	°C
Maximum operating non-condensing humidity	70	%

Degassing (per channel)

Parameter	Value	Unit
Degassing	Permeation through solid membrane	
Membrane material	Vinyl methyl silicone (VMQ)	
Degassing fluid volume	2	mL
Flow rate for 50% degassing efficiency (Water 25°C, Vacuum level 16.7kPa)	20 333	mL/min µL/sec
Recommended minimum degassing flow rate	5 83	mL/min µL/sec
Recommended maximum degassing flow rate	20 333	mL/min uL/sec
Recommended degassing vacuum	16.7	kPa
Vacuum volume	31	cm3
Pumpdown period to achieve 16.7 kPa vacuum using an IDEX Health & Science double stage vacuum pump	< 10	sec
Vacuum connection	4.7 mm OD	
Vacuum connection location	At end of tank	
Vacuum connection material	Polyphenylene sulfide (PPS)	
Recommended vacuum connection	Low gas permeability 3 mm ID elastomeric tube	

Regulatory

Parameter	Value
REACH	Yes
RoHS	Yes

Fluidic (per channel)

The following table contains nominal values for reference:

Parameter	Value	Unit
Internal fluid pathway	Parallel circular tubes	
Inline fluid volume	2.1	mL
Pressure drop (Water 25°C, Maximum operating flow rate)	< 20	kPa
Fluid contact materials	Ethylene propylene diene monomer (EPDM) Polyphenylene sulfide (PPS) Vinyl methyl silicone (VMQ)	
Fluid connection	Female flat-bottom 1/4-28UNF	
Fluid connection location	On top lid	
Fluid connection material	Polyphenylene sulfide (PPS)	
Fluid connection center-to-center distance	23.8	mm
Fluid connection edge-to-edge distance	14.3	mm
Recommended fluid connection	Male 1/4-28UNF nut preferably with a softer ETFE ferrule	

Mechanical

Parameter	Value	Unit
Housing material	Polyphenylene sulfide (PPS)	
Outer dimensions	75.7 x 49 x 22	mm
Mounting	With 6-32UNC screw	

Environmental

Parameter	Value	Unit
Operating temperature	5 to 40	°C
Operating non-condensing humidity	10 to 70	%
Storage temperature (The membrane must be dry inside.)	-10 to 60	°C
Storage non-condensing humidity	10 to 70	%
Vibration	Liquid delivery pumps with high vibration, such as diaphragm pumps, are not recommended	



Phone: 212-204-0075
Email: info@pyvot.tech
Web: www.pyvot.tech



**BIOTECH
FLUIDICS**

Typical Performance (per channel)





Figure 1: The plot shows the reference degassing efficiency of degassing chamber DGC-020 tested under various flow rates. The degassing efficiency represents the percentage of removed dissolved oxygen with respect to its saturation concentration. The curve shows a typical minimum performance profile gained at 25 °C using deionized water and a vacuum level of 16.7 kPa. The chamber was operated with negligible back pressure connected to the discharge side of a continuous pump. Oxygen content was monitored using a dissolved oxygen sensor while the flow rate was measured using a mass flow meter.


Integration (liquid integration shown for one channel)

The pump should ideally pull fluid through a degassing chamber to prevent from being pressurized above the maximum operating pressure of 100kPaG. It is not recommended to install the degassing chamber on the dispense side of the pump. Leave the cap on the unused vacuum port. Always place a liquid trap between the degassing chamber and a vacuum pump for longer continuous operation in order to avoid spillovers into the vacuum line.




Mounting Orientation

Condensation will reach an equilibrium in the chamber such that a small amount of water will be present at the bottom of the vacuum space. We recommend mounting the DGC-020 either in an upright position with fluid connections pointing upward and vacuum connections pointing downward or horizontally with fluid and vacuum connections pointing sideways with connecting elastomeric tube to the lower vacuum connection.



Dimensions



Warranty

Seller warrants to a buyer that each product will be free of defects in workmanship and material for the period of 1 year. The warranty period for all products commences on the date the product is deposited by the seller with the carrier for shipment. For complete warranty details refer to IDEX Health & Science LLC terms and conditions of sale which can be found at <https://www.idex-hs.com/about-us/legal-notice/terms-conditions-of-sale/>.



Phone: 212-204-0075
Email: info@pyvot.tech
Web: www.pyvot.tech



**BIOTECH
FLUIDICS**