




ヒノキチオール分離

Hinokitiol

Sunniest C18 5 μ m, 150 x 4.6 mm i.d. Metal free column



1. Hinokitiol



Column: Sunniest C18 5 μ m, 150 x 4.6 mm i.d. Metal free column (catalog number: EB347MH)

Mobile phase: Methanol:0.2 mM EDTA + 5 mM citrate buffer (pH4.0) = 40:60

Flow rate: 1.0 mL/min

Temperature: 40 $^{\circ}$ C

Detection: UV@240 nm

Sample: 1 = Hinokitiol